



Inactivation & Transport Fluid Kit

Cat. No: 23S-04-01

Ingredients: Tris-HCL , EDTA, Protein & Lipid spoiling

Physical Properties: Liquid, colorless

Description: The product is used for the transfer of enveloped and unenveloped viruses in a inactive condition for detection and research purposes.

Product Specification: The product is used in molecular virology, molecular biology, molecular genetic and chemistry fields. It is suitable for RNA/DNA isolation procedures and it does not contain DNase, RNase, exonucleases and endonucleases.

Only for laboratory use.

Procedure:

1. Immerse nasal/ nasopharyngeal swab sample directly in the inactivation solution and rotate for 3-5 times.
2. Discard the swab as medical waste.
3. Incubate inactivation solution at room temperature for 10 minutes.
4. Virally inactivated sample can be used for DNA/RNA isolation and all following molecular studies including Real Time PCR.

Cautions:

- In case of contact, it may irritate skin.
- Do not use without gloves.
- In case of contact, immediately wash skin with copious amounts of water.
- Do not use expired products.
- Shelf-life of the product is 18 months.
- Please contact us for your questions, tech@snp.com.tr

Storage: Room temperature

Date of issue: 01.04.2020

SNP Biotechnology R&D Ltd.

Hacettepe Technopolis 1. R&D Building Office No: 10 ANKARA / TURKEY

Tel.: +90.312.2992328 / Fax: +90.312.2992348 / www.snp.com.tr / info@snp.com.tr